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**GB - INSTRUCTIONS FOR USE**  
**EVOLUTIS instruments**  
**VERSION 01-2021**

For the attention of the surgeon and operation theatre staff

**Warnings**

EVOLUTIS instruments are to be exclusively used for the fitting or explantation of EVOLUTIS prosthesis. They are composed of: cannulated (or not) drills / reamers, tray, scalers, screwdrivers, rasps, burs, trial prosthesis, cutting guide, handles, extractors, impactors, grippers, guide pins for knee, hip, shoulder, extremities / traumatology surgery.

It is important that the surgeon has mastered the surgical technique and that they have fully read and understood these instructions and the operating technique manual specific to the selected joint implant. The operating technique manual is provided together with the material supplied to the establishment. Alternatively, it can be consulted online or downloaded from [www.evolutisfrance.com](http://www.evolutisfrance.com) upon the surgeon's request to Evolutis to obtain a personal access code for the download website. The instruments are not for single use. There is no known adverse effect when the instruments are used according to these instructions and the operating technique.

Only use cleaning agents and/or disinfectants specially formulated for cleaning or disinfecting reusable instruments.

The instructions provided by the detergent manufacturer regarding concentrations and temperatures must be followed.

**Limitations of Treatment**

These instruments have been designed for the fitting or explantation of an implant manufactured by EVOLUTIS; see the indications and contraindications of the latter.

**1- Initial treatment at the point of use**

Some instruments must be used with a motor, in this case the connection system must be a chuck or an AO tip or a Harris tip as appropriate. It is the responsibility of the operating theatre personnel to check, before the operation, the proper coupling of the instruments with the operating room equipment available. Some instruments are used to take measurements (in millimetres or degrees); the operator must ensure that the measurement is legible. Any damage to the equipment may result in malfunction.

The measuring accuracy for the lengths is  $\pm 1$  mm, and of  $\pm 2^\circ$  for the angles unless otherwise stated on the instruments. The instruments should only be used for the function assigned to them in the operating techniques.

In the event of a change in the performance of the instrument, it is imperative to inform the manufacturer and quarantine the instrument.

At the end of the operation, the surgeon should pay special attention to ensure that an instrument has not been left in situ.

Instruments must not be modified.

## **2- Preparation before cleaning**

Special attention should be given to instruments with several parts and/or with holes. Disassemble these if necessary.

## **3- Cleaning**

After any surgical procedure, the instruments should be cleaned as soon as possible (within 30 mins). Soak soiled instruments for at least 10 minutes in a mildly alkaline enzymatic detergent. Use validated cleaning procedures and rinse thoroughly (Example: products - DR WEIGERT). Do not use products such as soda.

## **4- Disinfection**

Soak soiled instruments for at least 10 minutes in a mildly alkaline enzymatic detergent (Example: products - DR WEIGERT). Use validated disinfection procedures and rinse thoroughly. Do not use products such as soda. In case of contact with patients at risk (NCTA) the instruments must be treated according to the procedures in force in the establishment.

## **5- Drying time**

30 minutes.

## **6- Maintenance, inspections and tests**

Instrument experience normal wear and tear during operations; when operating theatre personnel feel that an instrument is no longer usable, it must be decontaminated, cleaned and returned to EVOLUTIS.

For articulated instruments (e.g. pliers) the use of a lubricant (example: products - DR WEIGERT) validated for steam sterilisation processes can be used.

## **7- Sterilisation**

The instruments provided by Evolutis are non-sterile.

Before use, surgical personnel must ensure that the instruments have been sterilised. Instruments must be sterilised by trained and skilled personnel. Instrument boxes are not a sterile barrier: sterilisation packaging must be used to maintain the sterile conditions.

EVOLUTIS instruments have been validated for the following steam autoclave sterilisation cycle :

Cycle Method	Temperature	Duration
Wet heat sterilisation	132°C (269.5°F)	4 mn

It is essential that the sterilisation method is validated for the particular type of sterilisation equipment and for each load/product configuration according to ISO 15883.

### **8- Packaging**

Clean instruments must be packaged appropriately. Special attention must be paid to cutting instruments, which must be protected.

### **9- Storage**

Instruments must be stored in a clean, dry room at a temperature of 5 to 40°C (41 to 104°F).

### **10- Additional information**

In the event of a special or unusual risk, the instruments must be disposed of according to the legislation in force.

The various materials used for the instruments comply with the standards in force.

Metallic materials are radiopaque, the plastic materials are radiotransparent.

In the event of a severe incident, the user must notify the manufacturer [qualite@evolutis42.com](mailto:qualite@evolutis42.com) and the competent authority.

### **11- Manufacturer contact information**

For further information please contact EVOLUTIS (Tel 33 (0) 44 77 60 79 99) or its representative.

### **12- Composition**

- Stainless steel compliant with NF EN 10088-3, NFS 94-090, NF EN 10088-2 or NF EN 10028-7 or ASTM F899, NF EN 10216-5, ASTM A213 or NF EN 10151
- Stainless steel compliant with 5832-1 or 5832-9
- Cobalt chrome compliant with 5832-12 or 5832-4
- Titanium compliant with 5832-3 or 5832-2
- Polymers

*The instructions provided above have been validated by the manufacturer of the medical device for preparing a medical device for reuse. It remains the responsibility of the operator to ensure that the treatment, as it is actually applied by implementing the equipment, materials and personnel of the treatment facility, achieves the desired result. This requires verification and/or validation and routine inspection of the process.*

*These instructions have been drawn up in accordance with Annex B of ISO 17664-2017*

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