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For the attention of the surgeon and operation theatre staff

1 - GENERAL

The operating surgeon needs to be trained and familiar with the surgical technique, and must have read and understood the IFU (Instructions For Use) leaflet and the surgical technique manual provided by EVOLUTIS with the shipment of the material to the hospital, or consulted on-line or downloaded from the www.evolutisfrance.com site once the surgeon has requested and received a personal login code.

Equally to the surgical procedure and to the positioning of the implants, the selection of the patient is important.

An overweight or unsuitable functional requirement can create exceptional constraints and limit the lifespan of the arthroplasty.

The warnings and the instructions must be strictly observed. Precautions and instructions should be rigorously followed.

2 - INDICATIONS - CONTRA-INDICATIONS – COMPLICATIONS

Total knee arthroplasties are indicated for the treatment of symptomatic pain and/or functional problems of the knee in patients whose skeleton is mature and only when pain killer medication and correctly followed conservative treatment have failed. For the patient, his anatomy and the structure of his articulation will need to be adapted to receive the selected implant(s).

2-1-Indications

The indications for total knee arthroplasty are:

- Degenerative non inflammatory knee disease (gonarthritis).
- Inflammatory knee disease (rheumatoid arthritis, post traumatic arthritis).
- Metabolic knee disease (chondrocalcinosis).
- Post Traumatic degenerative arthritis.
- Functional repair of a recent trauma (fracture, dislocation)
- Revision after failure of a conservative surgery, a ligamentoplasty, a menisectomy, an osteosynthesis, or a partial or total arthroplasty.
- Periarticular or bone tumor.

2-2-Limits of indication

In primary surgery of the knee joint, and even more in revision or in tumoral surgery, an important deviation of the anatomic axis of the lower limb (genu-varum or genu-valgum superior to 15° of deviation), or the quality of the bone stock and the bone defects due to the explantation of any previously implanted material, can compromise the primary fixation

of the implantable device and thus limit its indications. Depending on the degree of axial deviation, or on the location and the extension of the bone defect, a knee implant with complementary fixation means or bone loss substituting elements such as intramedullary keels or wedges and augments, or with a more constraint design can be considered.

2-3-Contra-indications

The arthroplasty of the knee can be contra-indicated in cases of: local or systemic infection, mental deficiency, neuromuscular afflictions, vascular or neurological affections, patients subject to alcohol or psychotropic problems, excessive medication, excess functional situations (such as carrying out a sport where the risk of falling is prevalent or a situation in which the forces on the implant would exceed material resistance), overweight, insufficient bone stock or weak demineralized bone impeding a good fixation, or severe extra articular deformation.

Implantation of prosthesis is not recommended in pregnant women

2-4-Complications

Implanting a total knee prosthesis can lead to undesirable effects such as haematoma, thrombosis, pulmonary embolism, cardio-vascular disturbances, severe nervous, tendinous or venous risk factors, excessive peri-prosthetic ossification, allergic reaction to the material of the implant, tissue reaction to metal (metallosis) or to the polyethylene particles (granulomas), pain, bone fracture, breakage of the implant, bone wear, prosthetic wear, limb length discrepancy, implant loosening, dislocation, or infection.

3 - PRECAUTIONS

3-1-Pre operative phase

Before the surgery, the surgeon should discuss with the patient and before surgery, all physical and mental limitations particular to the patient and products for all aspects of the surgery at the pre, intra, and post-operative steps, including all the possible complications. This discussion should include limitations of the procedure and the mechanical limitations of the material(s) selected for implantation. Factors that could impose limits on the performance and stability of the implant(s) (ex: activity level and patient weight) should be described to the patient so as to maximize the potential for long-term freedom from complications.

The necessity to follow the surgeon's instructions regarding these factors following the operation must be clearly understood by the patient. The patient must be instructed that any change in the performance of the implant(s) needs to be immediately reported to the surgeon for proper action. The surgeon must make sure that the instructions have been clearly understood by the patient.

An adequate inventory of sterile implants sizes and the dedicated instrumentation set should be on hand at the time of surgery.

3-2-Manipulation

For sterile implants the integrity of the packaging must be verified as it is a guarantee of sterility, as well as end of shelf life date. If the packaging is damaged do not use, do not re-sterilize. Implants which have been unpacked and not implanted should not be re-sterilized. If implants are delivered non sterile they must be sterilized by recognized methods. Implants must be manipulated with caution in order to avoid damage or scratching. They should be conserved in their original sealed packaging between 5 and 40°C. If protective components are delivered with the implant, they must not be removed until implantation of the implant. Implants are designed for single use and must never be re-implanted.

3-3-Operative (Surgical) technique

The surgeon must be fully conversant with such surgical techniques. Careful pre-operative planning on the basis of radiographic findings should be carried out routinely. Radiographic templates for each implant in a variety of magnification ratio are available from EVOLUTIS.

3-4-Fixation of implants

Coated implants must not be implanted with orthopaedic bone cement. When cement is used, it is essential that the cement manufacturer's "Instructions for Use" be followed carefully. Handling of the bone cement can influence the efficacy of the implant fixation. Care must be taken that the correct and appropriate size implant is used in conjunction with the correct instrumentation and the trial components where these are provided. Definitive implants and trial components manufactured by Evolutis must not be used in conjunction with those of another manufacturer as components may not be compatible. The compatibility rule of the implants is described in the surgical technique of the implants through indications of sizes and match of the articular surfaces according to the type of postero-stabilization and to the type of tibial baseplate used (fixed or mobile).

It is important to identify clearly the implant before selection and opening by a complete reading of its labeling. The labeling includes indications of type of postero-stabilization, mode of fixation, side of implantation, type of tibial baseplate, size, and thickness.

3-5-Post operative follow up

The patient should be released from the hospital with complete instructions and warnings, preferably in writing, regarding any further exercises and therapies, and any limitation on his activities.

A continuing program of follow-up controls is essential and should be strictly enforced. In every case accepted practices should be followed in post-operative care of total knee arthroplasty.

In the event of a serious incident, the user must notify the manufacturer and the competent authority.

3-6-Removal of implant(s)

In case of removal of the implant, the surgeon must use accepted practices with the right instrumentation.

4 - WARNINGS

4-1-Re-use of an implant

The implants are intended for single use only; re-use of any implant is strictly prohibited. Should the implant be re-used its fixation and wear characteristics will be altered and the patient security put at risk.

In case of explantation of an implant, this implant should be eliminated in compliance to legislation in force.

In case of complete unpacking of an implant by mistake, this implant should be eliminated in compliance to legislation in force.

4-2-Imaging exam of a patient with an Evolutis prosthesis

Non-clinical testing has demonstrated that this prosthesis is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3T or less
- Maximum spatial gradient magnetic field of 19 T/m
- Maximum MR System reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the product is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

The appearance of an artifact caused by this device remains possible.

That's why the patient with a prosthesis must inform the healthcare staff before any MRI examination of the presence and site of implantation of his prosthesis.

5 - MATERIALS

The different materials used in knee implants are:

Cobalt chrome: ISO standard 5832-4 or 5832-12

Polyethylene PE-UHMW: ISO standard 5834-1 and 2

Stainless steel: ISO standard 5832-1

Titanium Ti6Al4V: ISO 5832-3

Pure titanium (coating): ISO standard 13179-1

Calcium Hydroxyapatite (coating): ISO standard 13779-2, 3 and 4

The material used for each individual device is specified on the labeling of each device.

6 - IMPLANT CARD

6-1-To the attention of medical staff

An implant card is provided in the prosthesis package. You should write the patient's name and the date of the operation and attach the patient label

to the space provided. This card must be given to the patient along with the information specified below.

6-2-For the attention of the patient

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during an examination such as MRI, CT scan or ultrasound.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

7 - FURTHER INFORMATION

For any further information please contact Evolutis or its local representative.

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