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VIS D'INTERFERENCE BIORESORBABLE

**GB** - INSTRUCTIONS FOR USE

ABSORBABLE INTERFERENCE SCREW

**ES** - MANUAL DE INSTRUCCIONES

TORNILLO DE INTERFERENCIA BIORREABSORBIBLE

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VITI AD INTERFERENZA

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**GB - INSTRUCTIONS FOR USE**  
**ABSORBABLE INTERFERENCE SCREW (FIXION)**  
**AND STAPLES**  
**VERSION 01-2021**

For the attention of the surgeon and operation theatre staff

## **1 - GENERAL**

It is of utmost importance that the surgeon be well informed on ligamentoplasty techniques and that he has read and paid particular attention to this instruction for use (IFU) leaflet. The surgical technique is delivered to the hospital together with the implants and the dedicated instrumentation set, or can be viewed directly, or downloaded under pdf format, on the [www.evolutisfrance.com](http://www.evolutisfrance.com) website after the surgeon has been awarded a personal access login.

In the same way that the implantation and positioning of the implant is important so is the patient selection criterion. Warnings and instructions must be rigorously followed.

## **2 - INDICATIONS - CONTRA-INDICATIONS – COMPLICATIONS**

### **2-1-Indications**

The essential indication of the interference screws is the fixation of the tendinous transplant in the surgical repair of the rupture of the ACL (anterior cruciate ligament) of the knee in cases of persistent instability of the knee (Lachmann's positive test, rotatory click) during intense occasional or high-level sports or professional practices, and when a past medical history (PMH) of meniscal lesion has created the knee instability.

The knee ACL ligamentoplasty repair can be indicated in cases of chronic instability in daily activities for non-sporting patients, or in cases of persistence of the instability in a conservative therapeutic initiative and only after a well driven physiotherapy has failed.

### **2-2-Contra-indications**

The young patient not having reached his skeleton maturity is an absolute contra-indication to the knee ACL ligamentoplasty repair.

The relative contra-indications to the knee ACL ligamentoplasty repair are: patient older than 50 years, patient not practicing intense occasional or high-level sports or professional practices, and clinical symptoms of sporadic instability not associated with true functional limitation.

The presence of an infection, an on-going medical treatment, a tumor, a mental deficit, a muscular disease, or a nervous or vascular infringement must be taken into account in the decision and in the timetable of management of the patient.

Implantation of prosthesis is not recommended in pregnant women

### **2-3-Complications**

As for any surgical operation, the ACL surgical repair can lead to undesir-

able effects such as failure of the ligamentary repair, hematoma, inflammatory reaction, or infection.

### **3 - PRECAUTIONS**

#### **3-1-Pre operative phase**

Before the surgery, the surgeon must question the patient about any difficulties either physical or mental the patient might suffer from, and must clearly inform the patient of the implant that will be used and of the different aspects of the surgical procedure. The information must emphasize on the limits of the surgery, and mention the mechanical limits of the implant as well as the material it is made of. It is also advisable to inform the patient of the factors which could limit the clinical outcomes and which could compromise the stability of the implant (e.g. failing to comply with the physiotherapy protocol, excessive and early solicitation of the repair...), in order to avoid any unnecessary complication. The surgeon should make sure that the recommendations and warnings regarding the post-operative behavior of the patient have been fully understood and accepted by the patient.

Before proceeding to the surgery, the surgeon or the theater staff must have checked that a sufficient number of implants of all available sizes are available.

#### **3-2-Manipulation**

Check that the temperature control indicator on the package of the implant has not changed color from BLUE to BLACK (BLUE: OK for use, BLACK: do not use). If the "marker dot" has turned BLACK do not use the screws.

Implants should be stored in their sealed original and intact packaging, and preserved at a temperature between 5°C and 30°C (screw) / 5°C and 40°C (staples). The integrity of the packaging should be checked to ensure that the sterility of the content is not compromised. The sterility date should also be controlled before use. If the packaging is damaged do not use nor re-sterilize. Extreme care must be taken while handling the implant so that the implant is not scratched or damaged in any way.

The implant should be implanted only with use of the instrument set delivered by Evolutis. The implants are made for single use only and should never be re-implanted.

#### **3-3-Operative (Surgical) technique**

The surgeon must use and follow documented and approved surgical techniques. A precise surgical planning procedure must be carried out in all cases based on clinical and medical imagery evaluations.

#### **3-4-Fixation of implants**

Care must be taken that the correct and appropriate size of implant is used in conjunction with the proper instrumentation. The implants supplied under an Evolutis brand must not be used with any instrument supplied by another manufacturer as these may not be compatible.

FIXION® interference screws can be used alone (2 per surgery) or in combination with other fixation devices (staples).

### **3-5-Post operative follow up**

The patient should be released from the hospital with complete recommendations and warnings, preferably in writing, regarding further exercises and therapies, and any limitations on his activities. A continuing program of follow up controls is essential and should be strictly enforced. In any case accepted practices should be followed in post-operative care.

In the event of a serious incident, the user must notify the manufacturer and the competent authority.

### **3-6-Removal of implant(s) and reesterilization**

In case of removal of the implant, the surgeon must follow documented and approved surgical techniques with the appropriate instrumentation.

## **4 - WARNINGS**

### **4-1-Re-use of an implant**

Implants are single use and their re-use is forbidden. In case of re-use, the implant stability and its wear can be altered and the patient security compromised. If an implant is unwrapped by mistake or removed, they must be disposed of according to the legislation.

### **4-2-Imaging exam of a patient with an Evolutis prosthesis**

Non-clinical testing has demonstrated that this prosthesis is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3T or less
- Maximum spatial gradient magnetic field of 19 T/m
- Maximum MR System reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the product is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

The appearance of an artifact caused by this device remains possible.

That's why the patient with a prosthesis must inform the healthcare staff before any MRI examination of the presence and site of implantation of his prosthesis.

## **5 - MATERIALS**

The interference screws are made of the following materials:

Bio-absorbable PDLA (D Lactide and L Lactide Poly Lactic Acid Co-polymer)

Staple material: cobalt-chromium alloy in compliance with ISO 5832-4.

## **6 - IMPLANT CARD**

### **6-1-For the attention of the medical personnel**

An implant card is provided in the prosthesis package. You should write

the patient's name and the date of the operation and attach the patient label to the space provided. This card must be given to the patient along with the information specified below.

### **6-2-For the attention of the patient**

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during an examination such as MRI, CT scan or ultrasound.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

### **7 - FURTHER INFORMATION**

For any further information please contact Evolutis or its local representative.

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