## 1 – GENERALITÀ

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## 1 – GENERAL

This instruction for use is intended to apply to the TMC HORUS and ISIS product line, including following components:

- Metacarpal stems, coated with titanium plasma spray, with 5 sizes
- Metacarpal modular head, with offset neck 15°, in size M and L
- Metacarpal modular head, with straight neck, in size M and L
- Threaded Trapezium Cups, cementless, with titanium plasma spray, in Size M and L
- Pexel-E liners for threaded trapezium cups, in size M and L
- Pexel-E cemented trapezium cups, in size M and L

The TMC joint arthroplasty requires qualified surgeons for such a surgery, being aware and having understood the following warnings and being thoroughly familiar with the TMC operating technique.

It is mandatory to use the TMC instrumentation provided by EVOLUTIS.

Patient selection is as important as the surgical procedure. Following warnings and instructions must be scrupulously observed.

## 2 - INDICATIONS - CONTRAINDICATIONS - COMPLICATIONS

### 2.1 - Indications

The general indications for the use of a TMC joint prosthesis are the following:

- The presence of pain is the fundamental element for deciding on surgery. Pain can be the result of traumatic or degenerative pathologies.
- Post-traumatic sequelae
- Failure of previous joint surgery

The TMC joint prosthesis is indicated for:

- A primary 3rd or 4th stage arthritis, or post-traumatic arthritis of the basal thumb joint (rhizarthrosis)
- patient over the age of 60, without any heavy manual activity, nor in professional practices, in leisure time or in sport.
- Minimum thickness of the trapezium bone > 8 mm

### 2.2 Intended use

The TMC trapeziometacarpal implants are intended to replace the natural joint between the trapezium and the first metacarpal of the thumb.

The HORUS and ISIS joint prosthesis are indicated for replacement in case of advanced destruction of the TMC joint.

For this purpose, the HORUS and ISIS joint prosthesis will be implanted in the thumb and trapezium and will come into contact with bones and tissues permanently.

The HORUS and ISIS prosthesis are composed of a metacarpal implant, a trapezium cup and in between, a modular articular head that allows to articulate the two components and adjust the stability of the joint.

### 2.3 Intended User

It will be compulsorily implanted by validated and trained orthopedic surgeons. The HORUS and ISIS TMC surgical technique shall be applied

### 2.4 - Contraindications

Arthroplasty can be contraindicated in the following cases:

- Acute or chronic, local or systemic infections.
- Serious muscular, neurological or vascular disabilities affecting the extremity in question (great caution is recommended in case of inflammatory rheumatism)
- Bone degeneration or poor bone quality that can affect the stability of the implant
- Any concurrent conditions that may affect the function of the system
- Status of the ligament system unsuitable for the type of prosthesis chosen.
- A shallow trapezium (< 8mm) is a severe contraindication to the trapezio-metacarpal prosthesis.

### 2.5 Formal contraindication

- Inappropriate ligament function
- Poor bone quality or bone stock that can affect implant stability.
- Trapezium must be deep enough to allow the placement of cup in good conditions of stability and without degrading the cortical walls

### 2.6 – Possible complications

The implantation of a prosthesis can generate undesirable effects such as hematoma, thrombosis, pulmonary embolism, cardio-vascular disorders, serious nervous, tendon or venous risk factors, excessive periprosthetic ossification, allergic reaction to the implant material, tissue reaction to metallic particles (metallosis), pain, bone fracture, breakage of the prosthesis, bone wear, wear of the prosthesis, loosening, dislocation or infection of the prosthesis.

## 3-PRECAUTIONS

### 3-1 - Preoperative phase

Before surgery, the surgeon must investigate about patient physical and mental state and his difficulties and inform him about the product and the different aspects of the operation. During the visit, the surgeon must explain the limits of the intervention and the mechanical or material limitations of the selected prostheses.

Factors that can limit the performance of the prosthesis and alter its stability (e.g. functional overload) must be explained in order to reduce and delay the onset of complications related to them. The need to follow the surgeon's instructions and prescription should be explained to the patient. The appropriate prostheses must be available at the time of the operation.

### 3-2 - Handling

As implants are gamma sterilized, the integrity of package must be checked before use, which guarantees the sterility of content. The expiry date of each individual package must be checked before opening.

If the package is damaged, do not use or steam re-sterilize; Return the implant in the defective package to EVOLUTIS or its representative.

Prostheses that have been opened and not implanted will not be re-sterilized. The prostheses are designed for single use and must not be re-implanted.

Prostheses should be handled with care to avoid damage or scratches. They must be kept in their original sealed packaging between 5 and 30°C.



Instruction for opening the packaging without risk :

External Casing


External Casing Opening : remove each corner first, one by one

Internal Casing Display

Internal Casing Opening

### 3-3- Operating technique (surgical)

Careful preoperative planning, based on radiological examination and cast to determine the size of the prosthesis to be implanted, must be carried out systematically.

### 3-4 - Fixing the prosthesis

Coated prosthesis shall not be implanted with cement. When using cement for the cup, it is essential to follow the cement manufacturer's instructions and advice. The way of cement is prepared can influence the stability of trapezium cup.

An appropriately sized prosthesis with appropriate provided instrumentation and trials shall be used.

The HORUS and ISIS components and the trials supplied by EVOLUTIS shall not be used with components supplied by another manufacturer as they are not compatible. The compatibility of the final prostheses is mentioned in the surgical technique with information about diameter, size and length.

### 3-5 - Post-operative phase

Towards the end of the hospital stay, patients shall receive written advice and warnings on complementary treatments, ie. physiotherapy and activity restrictions. Regular monitoring with periodic checks is essential and shall be strictly applied.

### 3-6 - Removal of prosthesis

In the event of prosthesis removal need, the surgeon must use proven surgical techniques and adequate Evolutis instrumentation.

## 4 – Warnings

### 4.1 - Reuse of prosthesis

The reuse of HORUS and ISIS components is prohibited. In case of reuse, the stability of the prosthesis and its wear could be altered, and the safety of the patient compromised. Prostheses mistakenly opened (no longer sterile) or explanted must be disposed of, according to current legislation.

### 4.2 - Magnetic resonance environment warning

The HORUS and ISIS prosthesis have not been tested in an electromagnetic field (magnetic resonance imaging) due to the fact that the materials used for these implants are little or no ferromagnetic and without electrical activity.

More generally, the presence of an Evolutis metacarpal prosthesis is not a contraindication for magnetic resonance imaging (Max 3 Tesla). However, the following potential risks exist:

- Artifact risk: Magnetic resonance imaging can produce artifacts that alter the interpretation of the imaging exam on the periphery of the prosthesis. However, outside the immediate vicinity of the prosthesis, the anatomical structures become perfectly analyzable again.
- Risk of overheating: MRI can cause local heating.
- Risk of displacement or twisting: Magnetic resonance imaging can cause displacement or twisting.

The patient with a prosthesis must inform healthcare professionals prior to any MRI examination of the presence and implantation site of the prosthesis.

### 4.3 Logistic advices

Warning: the packaging allows air and land transport

Before surgery, the completion of implant set and instrument set must be guaranteed, according to the checklist of components.

Restriction on use : No mix is allowed between TMC ISIS or HORUS component and TMC components from other manufacturer

## 5- MATERIALS

The different materials used with HORUS and ISIS prosthesis are: Stainless steel: ISO 5832-9 and ISO 5832-1 standards

Cobalt chrome: ISO 5832-12 standard

Titanium alloy: ISO 5832-2 and ISO 5832-3

Polyethylene: ISO 5834-1 and -2

For coatings: Pure titanium ISO 13179-1 e ASTM F1580 and titanium nitride

6- For further information, complaint, feedback and local correspondence for vigilance, please contact :

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 Manufacturer	 Catalog number	 Production Batch N°	 Expiry date	 Sterile	 Store between 5°C and 30°C
 Do not reuse	 Caution, consult Operative technique	 Consult instructions for use	 Do not use if package is damaged	 Keep dry	 Keep away from sunlight

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**EVOLUTIS SAS**