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**GB - INSTRUCTIONS FOR USE**  
**OSTEOSYNTHESIS DEVICES FOR SHOULDER TRAUMA**  
**VERSION 2023-01**

For the attention of the surgeon and operation theatre staff

## **1 - GENERAL**

The operating surgeon needs to be trained and familiar with the surgical technique, and must have read and understood the IFU (Instructions For Use) leaflet and the surgical technique manual provided by EVOLUTIS with the shipment of the material to the hospital, or consulted on-line or downloaded from the [www.evolutisfrance.com](http://www.evolutisfrance.com) site once the surgeon has requested and received a personal login code.

Equally to the surgical procedure and to the positioning of the implants, the selection of the patient is important.

An unsuitable functional requirement can create exceptional constraints and limit the lifespan of the arthroplasty.

The warnings and the instructions must be strictly observed. Precautions and instructions should be rigorously followed. It is specified that all implants marketed by EVOLUTIS are supplied sterile (sterilisation by gamma radiation).

## **2 - INDICATIONS - CONTRA-INDICATIONS - COMPLICATIONS**

### **2-1 - Indications**

Shoulder osteosynthesis implantable devices are indicated for the treatment of acute and fresh trauma of the proximal humerus. The JUST UNIC device is indicated for Neer 3 and 4 comminuted fractures of the proximal humerus (3 or 4 parts) in osteoporotic bone for patients above 55 years old, and in fracture-dislocation cases with autologous bone-graft for patients under 55 years old.

The anatomy and structure of the patient's articulation will need to be adapted to receive the selected implantable device(s).

### **2-2 - Contra-indications**

Absolute contra indications of the JUST UNIC osteosynthesis are: shoulder trauma with extension to the humeral diaphysis, omarthrosis of the shoulder joint, deficient or unrepairable rotator cuff, revision of a failed shoulder arthroplasty, local or systemic infection, mental deficiency, neuromuscular disease, neurological or vascular problems, and patients addicted to alcohol or psychotropic drugs.

Relative contra indications of shoulder osteosynthesis with a JUST UNIC are: fracture at the young non-osteoporotic patient, excessive functional demands (sport with risk of falls or with excessive functional demands beyond the resistance limits of the device), overweight, insufficient bone stock or decalcification of the bone which would compromise the prosthetic fixation, severe humerus deformity, drug abuse, pre-existing periarticular oncological pathology.

Implantation of prosthesis is not recommended in pregnant women

### **2-3 - Complications**

Specific undesirable side effects may be associated with the JUST UNIC device: persistent pain, avascular necrosis/osteonecrosis, bone loss/resorption, mal-union/non-union of the fractured parts, peri-prosthetic ossification, allergic reaction to the material of the device, bone fracture, breakage or failure of the osteosynthesis and/or of the device, bone resorption, loosening of the device, dislocation.

General undesirable side effects may be associated with the surgical procedure: infection, hematoma, thrombosis, pulmonary embolism, cardiovascular disturbances, nerve stretching or injury resulting in alteration of sensitivity (muscle tingling, tenderness defects) or of motricity (reduction or loss of strength).

#### **2-4 - Benefits and performances**

The intended clinical benefits and performances of the JUST UNIC device are:

- An improvement of the success of the osteosynthesis associated with an improvement of the mobility of the shoulder,
- The healing of the tuberosities in frequency and position and the absence of secondary migration.

### **3 - PRECAUTIONS**

#### **3-1 - Pre operative phase**

Prior to the operation, the surgeon must interview the patient about his physical and mental status and difficulties and inform him about the implantable device and the different aspects of the operation. During the consultation the surgeon should explain the limits of the operation as well as the mechanical or the material limitations of the chosen device(s). Factors that could impose limits on the performance and stability of the implant(s) (ex: excessive functional expectation) should be described to the patient so as to maximize the potential for long-term freedom from complications, and the patient must be instructed that any change in the performance of the implant(s) needs to be immediately reported to the surgeon for proper action.

The necessity to follow the surgeon's instructions regarding these factors following the operation must be clearly understood by the patient. The surgeon must make sure that the instructions have been clearly understood by the patient.

#### **3-2 - Manipulation**

For sterile implants the integrity of the packaging must be verified as it is a guarantee of sterility, as well as end of shelf life date. If the packaging is damaged do not use, do not re-sterilize. Implants which have been unpacked and not implanted should not be re-sterilized.

Implants must be manipulated with caution in order to avoid damaging or scratching. They should be conserved in their original sealed packaging between 5 and 40°C. If protective components are delivered with the implant, they must not be removed until implantation of the implant. Implants are designed for single use and must never be re-implanted.

#### **3-3 - Operative (Surgical) technique**

The surgeon must use documented and validated surgical techniques. Careful pre-operative planning on x rays should be done systematically. Templates are available for the majority of the implantable devices at EVOLUTIS.

### **3-4 - Fixation of Devices**

Cementless devices must not be implanted with cement.

The appropriate device of the appropriate size should be used with the appropriate instrumentation and trials (if the latter are supplied).

### **3-5 - Post operative follow up**

At the end of their hospitalization, the patient must be provided with the precautions and prescriptions, preferably in writing, for the complementary treatments, physiotherapy exercises and restrictions on activity as required. Regular monitoring with periodic controls is essential and must be strictly applied. In all cases, the proven techniques for postoperative monitoring should be followed.

In the event of an MRI or CT scan, ask the patient to always disclose the fact that he/she has undergone an operation.

In the event of a severe incident, the user must notify the manufacturer and the competent authority

### **3-6 - Removal of the device(s)**

In case of device removal, the surgeon must use documented and validated surgical techniques and use the appropriate instrumentation

## **4 -WARNINGS**

### **4-1 - Re-use of an implant**

Implants are single use and their re-use is forbidden. In case of re-use, the implant stability and its wear can be altered and the patient security compromised. If an implant is unwrapped by mistake or removed, they must be disposed of according to the legislation.

### **4-2 - Imaging exam of a patient with an Evolutis prosthesis**

Non-clinical testing has demonstrated that this prosthesis is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3T or less
- Maximum spatial gradient magnetic field of 19 T/m
- Maximum MR System reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the product is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

The appearance of an artifact caused by this device remains possible.

That's why the patient with a prosthesis must inform the healthcare staff before any MRI examination of the presence and site of implantation of his prosthesis.

### **4-3 - Compatibility of implants and restrictions for use**

Implants and trials provided by Evolutis should not be used with components supplied by another manufacturer since they might not be compatible.

## 5 - MATERIALS

The different materials used in JUST UNIC devices are:

- Titanium alloy: ISO standard 5832-3 (stem, dock, staple, and screw)

ISO	C	N	Al	O	V	Fe	H	Other
5832-3	<0,08	<0,05	5.5 to 6.75	<0,13	3.5 to 4.5	≤0.30	≤0.015	Ti

- Cobalt-chromium: ISO standard 5832-12 (humeral head)

ISO	C	Cr	Ni	Mo	Mn	N
5832-12	≤0,14	26 to 30	≤1,0	5 to 5	≤1,0	<0,25

The material used for each individual device is specified on the packaging of each device.

## 6 - IMPLANT CARD

### 6-1 - For the attention of the medical personnel

An implant card is provided in the prosthesis package. You should write the patient's name and the date of the operation and attach the patient label to the space provided. This card must be given to the patient along with the information specified below.

### 6-2 - For the attention of the patient

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during an examination such as MRI, CT scan or ultrasound.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

## 7 - FURTHER INFORMATION

The summary of safety and clinical performance of device shall be made available to the public via Eudamed.

For any further information please contact Evolutis or its local representative.

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