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## **GB HIP : Stem / Head / accessories and screw**

*Amethyst / Cemtiv / Emergence / Evidence / Hactiv / Hagap / Helmed / Prius / Stemsys / X'Pert*

### **INSTRUCTIONS FOR USE** **VERSION 2022-05**

For the attention of the surgeon and operation theatre staff

#### **1 - GENERAL**

It is essential that the surgeon be fully conversant with the prescribed surgical technique and that due consideration be given to the information provided herewith. The specific surgical technique is delivered together with the implants and instruments to the care center before the surgery, or may be consulted online or downloaded from the [www.evolutisfrance.com](http://www.evolutisfrance.com) website. The surgeon will receive a personal login to access the website once he has registered on the website. Improper patient selection, with particular regard to patient weight and functional demands, as well as placement and positioning of the implant may result in unusual stress conditions and subsequent reduction in their functional life. Warnings and instructions should be strictly observed.

#### **2 - INDICATIONS - CONTRA-INDICATIONS - COMPLICATIONS**

Hemi and total hip arthroplasties are indicated for the treatment of symptomatic pain and/or functional problems of the hip in patients whose skeleton is mature and only when pain killer medicamentation and correctly followed conservative treatment have failed. For the patient, his anatomy and the structure of his articulation will need to be adapted to receive the selected implant(s).

##### **2-1-Indications**

Indications for hip arthroplasty are:

- Degenerative non inflammatory hip disease (coxarthrosis, arthritis of the hip).
- Inflammatory hip disease (rheumatoid arthritis, post traumatic arthritis).
- Metabolic hip disease (chondrocalcinosis).
- Post Traumatic degenerative arthritis.
- Avascular necrosis.
- Congenital dysplasia of the hip.
- Functional repair of a recent trauma (fracture, dislocation)
- Revision of a failed conservative surgery, osteosynthesis, partial or total arthroplasty.
- Tumoral surgery when affecting the hip joint.

The recommended indication for STEMSYS gritblasted are the displaced femoral neck fracture (Garden III or IV) of patients over 75 years of age.

Depending on the indication, their benefit/risk ratio is only evaluated in the short term.

Due to its shortening, the STEMSYS MI stem is not recommended for partial arthroplasty nor for revision arthroplasty.

The trochanteric grip plate is used in combination with the PRIUS stem and cortical screws Ø 4.5mm (see surgical technique) in case of:

- Fractures of the trochanteric region following the implantation of a femoral stem.

- Trochanterotomy for implantation or removal of a femoral stem (maintenance of the trochanter or femoral component).  
Treatment of trochanter pseudarthrosis after primary prosthesis or revision prosthesis.

In primary surgery of the hip joint, and even more in revision or tumoral surgery, the quality of the bone stock and the bone defects due to the ablation of any previously implanted material can compromise the primary fixation of the implantable device and thus limit its indications. Depending on the location and the extension of the bone defect, a longer cemented or cementless femoral component including a variety of complementary fixation means can be considered.

### **2-2-Contra-indications**

The arthroplasty of the hip can be contra-indicated in cases of: local or systemic infection, mental deficiency, neuromuscular disease, neurologic or vascular problems, patients addicted to alcohol or psychotropic drugs, drug abuse, excessive functional expectations (sport with risk of falls or with excessive functional expectations beyond the limits of the mechanical resistance of the prosthesis), overweight, insufficient bone stock or demineralisation of the bone which would compromise the prosthetic fixation, or severe extra articular deformity. Implantation of prosthesis is not recommended in pregnant women.

### **2-3-Complications**

The implantation of a hip prosthesis may entail the following complications: Haematoma, thrombosis, pulmonary embolism, cardiovascular disorders, lack of bone union, muscle, nerve, tendon or venous damage, peri-prosthetic ossifications, allergic reactions to prosthetic materials, tissue reaction to metal or plastic particles (metallosis-osteolysis), release of ions or metal debris, disassembly of prosthetic modules, displacement of the implant, disassembly of prosthetic modules, pain in the prosthesis, bone fractures, fracture of a prosthetic element, wear in the bones, wear of a prosthetic element, prosthetic noise (squeaking), uneven limb length, dislocation, loosening of the prosthesis, infection, acute renal failure, stroke.

## **3 - PRECAUTIONS**

### **3-1-Preoperative step**

The surgeon should discuss all physical and mental limitations particular to the patient and products for all aspects of the surgery at the pre, intra, and post-operative steps, including all the possible complications, with the patient before surgery. This discussion should include limitations of the procedure and the mechanical limitations of the material(s) selected for implantation. Factors that could impose limits on the performance and stability of the implant(s) (ex: activity level and patient weight) should be described to the patient so as to maximise the potential for long-term freedom from complications, and the patient must be instructed that any change in the performance of the implant(s) needs to be immediately reported to the surgeon for proper action.

The necessity to follow the surgeon's instructions regarding these factors following the operation must be clearly understood by the patient. The surgeon must make sure that the instructions have been clearly understood by the patient.

An adequate inventory of sterile implants sizes and the dedicated instrumentation set should be on hand at the time of surgery.

### **3-2-Handling of the implantable device(s)**

For implants that are provided sterile through Gamma irradiation, the integrity of the packaging should be checked to ensure that the sterility of the contents is not compromised. The expiry date of each individual packaging should also be controlled before opening. If the packaging is damaged do not use or re-sterilise, and return the implant together with its defective packaging to Evolutis or to its local representative. Unused components should not be re-sterilized.

Extreme care must be taken while handling the implant so that it is not scratched or damaged in any way.

Implants should be stored in their original sealed packaging and in a dry and tempered room (5-40°C). Where protective covers are provided these should be left in place until immediately prior to use.

Implants are designed for single use and must never be re-implanted.

### **3-3-Surgical technique**

The surgeon must be fully conversant with such surgical techniques.

Careful pre-operative planning on the basis of radiographic findings should be carried out routinely. Radiographic templates are available for the majority of implants.

### **3-4-Implant fixation**

Coated implants must not be implanted with cement.

When cement is used, it is essential that the cement manufacturer's "Instructions for Use" be followed carefully. Handling of the bone cement can influence the efficacy of the implant fixation.

Care must be taken that the correct and appropriate size implant is used in conjunction with the correct instrumentation and the trial components where these are provided.

The implantation of components assembled with a Morse taper connection require extra care from the surgeon especially if one of the components is made of ceramic material. The Morse taper on the femoral stem side must be thoroughly cleaned and dried before assembly of the components. If the component is made of ceramic, turn it clock-wise before impaction. Any impaction of a ceramic component must be made through a plastic instrument. Never impact a ceramic component directly with a metal tool.

### **3-5-Post operative follow up**

Towards the end of the hospital stay patients should receive written advice and warnings about complimentary treatment, physio-therapy exercises and activity restrictions

if necessary. In all cases validated post-operative follow up techniques must be used. The patient must be instructed to systematically inform that he or she has had an operation in case of subsequent MRI or scanner examination.

In the event of a serious incident, the user must notify the manufacturer and the competent authority.

### **3-6-Removal of implant(s)**

In case of removal of the implant, the surgeon must use accepted practices with the right instrumentation.

## **4 - WARNINGS**

### **4-1-Re-use of an implant**

Implants are single use and their re-use is forbidden. In case of re-use, the implant stability and its wear can be altered and the patient security compromised. If an implant is unwrapped by mistake or removed, they must be disposed of according

to the legislation.

#### **4-2-Imaging exam of a patient with an Evolutis prosthesis**

Non-clinical testing has demonstrated that this prosthesis is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3T or less
- Maximum spatial gradient magnetic field of 19 T/m
- Maximum MR System reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the product is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

The appearance of an artifact caused by this device remains possible.

That's why the patient with a prosthesis must inform the healthcare staff before any MRI examination of the presence and site of implantation of his prosthesis.

#### **4-3-Compatibility of implants and restrictions for use**

There may be restrictions on the use of certain implant sizes or models, particularly depending on their maximum mechanical strength. These restrictions shall be specified on the labels of the implants concerned.

The combination of a metal head with a conventional polyethylene insert is only recommended for patients over 70 years of age, or for patients under 70 years of age with a low life expectancy.

EVOLUTIS test implants and definitive implants must not be used with elements from any other manufacturer, as the components may not be compatible.

The compatibility of the implants is described in the surgical technique with indications for the Morse taper (combination of head/stem).

These characteristics are indicated on the implant labels.

Do not combine a metal head with a ceramic insert or a ceramic head with a metal insert.

Do not combine an extra-long neck head ( $\varnothing 28 + 7.0\text{mm}$  or  $\varnothing 32 + 8.0\text{mm}$  or  $\varnothing 36 + 8.0\text{mm}$ ) or a super-long neck head ( $\varnothing 28 + 10.5\text{mm}$ ) with a lateralised femoral stem.

Do not combine a skirted head ( $\varnothing 28 + 7.0\text{mm}$  or  $\varnothing 28 + 10.5\text{mm}$ ) with a dual mobility cup.

#### **4-4-Specific warnings in case of breakage of a ceramic component**

Whenever a ceramic component breaks – except if it happens in the intraoperative step– is deemed to be destructive for the metallic taper on which it is fixed. Consequently:

- In case of breakage of the ceramic femoral head when used in a ceramic-polyethylene articulation, the general guidance is to replace the femoral stem and the broken femoral head by a new femoral stem and a new ceramic head, except if – after a thorough inspection of the connection taper of the femoral stem – the taper is not damaged; in this case, the replacement of the femoral head alone may be considered. In any case, it is recommended to change the polyethylene insert.

- In case of breakage of one of the ceramic components when used in a ceramic-ceramic articulation, the general guidance is to replace the 2 ceramic articular components, the femoral stem and the acetabular cup. An exception can be made for the articular component which is not damaged if – after removal of the ceramic component and thorough inspection of the connection taper – the taper is not damaged: the replacement of the ceramic component alone can be considered.

- In case of breakage of the two ceramic components when used in a ceramic-ceramic articulation, the general guidance is to replace the 2 ceramic articular components, the femoral stem and the acetabular cup.

In all cases of breakage of a ceramic component, it is recommended to thoroughly clean the implantation site in order to remove any remaining ceramic fragment likely to induce third body complication.

## **5 - MATERIALS**

The different materials used in hip implants are:

- Stainless steel: ISO standard 5832-1 or 5832-9
- Titanium alloy: ISO standard 5832-2 or 5832-3
- Cobalt chrome: ISO standard 5832-4 or 5832-12
- Ceramic: ISO standard 6474
- Pure titanium (coating): ISO standard 13179-1
- Calcium Hydroxyapatite (coating): ISO standard 13779-2, 3 and 4

The material used for each individual device is specified on the packaging of each device.

## **6 - IMPLANT CARD**

### **6-1-For the attention of the medical personnel**

An implant card is provided in the prosthesis package. You should write the patient's name and the date of the operation and attach the patient label to the space provided. This card must be given to the patient along with the information specified below.

### **6-2-For the attention of the patient**

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during an examination such as MRI, CT scan or ultrasound.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

## **7 - FURTHER INFORMATION**

For any further information please contact Evolutis or its local representative.

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