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**GB - INSTRUCTIONS FOR USE**  
**ABSORBABLE INTERFERENCE SCREW (Fixion)**  
**VERSION 2023-10**

For the attention of the surgeon and operation theatre staff

## **1 - GENERAL**

It is of utmost importance that the surgeon be well informed on ligamentoplasty techniques and that he has read and paid particular attention to this instruction for use (IFU) leaflet. The surgical technique is delivered to the hospital together with the implants and the dedicated instrumentation set, or can be viewed directly, or downloaded under pdf format, on the [www.evolutisfrance.com](http://www.evolutisfrance.com) website after the surgeon has been awarded a personal access login.

In the same way that the implantation and positioning of the implant is important so is the patient selection criterion. Warnings and instructions must be rigorously followed. It is specified that all implants marketed by EVOLUTIS are supplied sterile (sterilisation by gamma radiation).

## **2 - INDICATIONS - CONTRA-INDICATIONS – COMPLICATIONS**

### **2-1-Indications**

The essential indication of the interference screws is the fixation of the tendinous transplant in the surgical repair of the rupture of the ACL (anterior cruciate ligament) of the knee in cases of persistent instability of the knee (Lachman's positive test, rotatory click) during intense occasional or high-level sports or professional practices, and when a past medical history (PMH) of meniscal lesion has created the knee instability.

The knee ACL ligamentoplasty repair can be indicated in cases of chronic instability in daily activities for non-sporting patients, or in cases of persistence of the instability in a conservative therapeutic initiative and only after a well driven physiotherapy has failed.

FIXION® interference screws can be used alone or in combination with other fixation devices (e.g. staples, cross-pin or endobuttons). One or two screws can be used during the procedure, the choice is left to the surgeon depending on the technique used by each operator.

### **2-2-Contra-indications**

The young patient not having reached his skeleton maturity is an absolute contra-indication to the knee ACL ligamentoplasty repair.

The relative contra-indications to the knee ACL ligamentoplasty repair are: patient older than 50 years, patient not practicing intense occasional or high-level sports or professional practices, and clinical symptoms of sporadic instability not associated with true functional limitation.

The presence of an infection, an on-going medical treatment, a tumor, a mental deficit, a muscular disease, or a nervous or vascular infringement

must be taken into account in the decision and in the timetable of management of the patient.

Implantation of implant is not recommended in pregnant women

### **2-3-Complications**

As for any surgical operation, the ACL surgical repair can lead to undesirable effects such as failure of the ligamentary repair, hematoma, inflammatory reaction, or infection.

### **2-4-Benefits and performances**

The intended clinical benefits and technical performances of the interference screws are:

- The ACL reconstruction using Kenneth-Jones technique or semitendinosus autograft and correction of Lachman and pivot-shift positive tests,
- The fixation of the interference screws during the first months post-surgery,
- The return to sportive activities after 6 months of traditional rehabilitation.

## **3 - PRECAUTIONS**

### **3-1-Pre operative phase**

Before the surgery, the surgeon must question the patient about any difficulties either physical or mental the patient might suffer from, and must clearly inform the patient of the implant that will be used and of the different aspects of the surgical procedure. The information must emphasize on the limits of the surgery, and mention the mechanical limits of the implant as well as the material it is made of. It is also advisable to inform the patient of the factors which could limit the clinical outcomes and which could compromise the stability of the implant (e.g. failing to comply with the physiotherapy protocol, excessive and early solicitation of the repair...), in order to avoid any unnecessary complication. The surgeon should make sure that the recommendations and warnings regarding the post-operative behavior of the patient have been fully understood and accepted by the patient.

Before proceeding to the surgery, the surgeon or the theater staff must have checked that a sufficient number of implants of all available sizes are available.

### **3-2-Manipulation**

Check that the temperature control indicator on the package of the implant has not changed color from WHITE to BLACK (WHITE: OK for use, BLACK: do not use). If the "marker dot" has turned BLACK do not use the screws.

Implants should be stored in their sealed original and intact packaging, and preserved at a temperature between 5°C and 40°C. The integrity of the packaging should be checked to ensure that the sterility of the content is not compromised. The sterility date should also be controlled before use. If the packaging is damaged do not use nor re-sterilize. Extreme care must

be taken while handling the implant so that the implant is not scratched or damaged in any way.

The implant should be implanted only with use of the instrument set delivered by EVOLUTIS. The implants are made for single use only and should never be re-implanted.

### **3-3-Operative (Surgical) technique**

The surgeon must use and follow documented and approved surgical techniques. A precise surgical planning procedure must be carried out in all cases based on clinical and medical imagery evaluations.

### **3-4-Fixation of implants**

Care must be taken that the correct and appropriate size of implant is used in conjunction with the proper instrumentation. The implants supplied under an EVOLUTIS brand must not be used with any instrument supplied by another manufacturer as these may not be compatible.

FIXION® interference screws can be used alone or in combination with other fixation devices (e.g. staples, cross-pin or endobuttons). One or two screws can be used during the procedure, the choice is left to the surgeon depending on the technique used by each operator.

### **3-5-Post operative follow up**

The patient should be released from the hospital with complete recommendations and warnings, preferably in writing, regarding further exercises and therapies, and any limitations on his activities. A continuing program of follow up controls is essential and should be strictly enforced. In any case accepted practices should be followed in post-operative care.

In the event of a serious incident, the user must notify the manufacturer and the competent authority.

### **3-6-Removal of implant(s) and resterilization**

In case of removal of the implant, the surgeon must follow documented and approved surgical techniques with the appropriate instrumentation.

## **4 - WARNINGS**

### **4-1-Re-use of an implant**

Implants are single use and their re-use is forbidden. In case of re-use, the implant stability and its wear can be altered and the patient security compromised. If an implant is unwrapped by mistake or removed, they must be disposed of according to the legislation.

## **5 - MATERIALS**

The interference screws are made of the following materials:

Bio-absorbable PDLA (D Lactide (30%) and L Lactide (70%) Poly Lactic Acid Co-polymer)

Material	L-Lactide polymer (%)	D,L-Lactide polymer (%)	Residue of Lactide monomer (%)	Residue of toluene solvent (ppm)	Tin (ppm)	Rate of sulphated ash (%)	Inherent viscosity (dl/g) before injection
<b>PDLLA</b>	67 - 73	27 - 33	< 0.1	< 890	< 60	< 0.1	5.7 – 6.5

The device:

- does not incorporate a substance which, if used separately, may be considered to be a medicinal product, including a medicinal product derived from human blood or human plasma.
- is not manufactured utilising tissues or cells of human or animal origin, or their derivatives.
- does not contain nanomaterials.

## 6 - IMPLANT CARD

### 6-1-For the attention of the medical personnel

An implant card is provided in the prosthesis package. You should write the patient's name and the date of the operation and attach the patient label to the space provided. This card must be given to the patient along with the information specified below.

### 6-2-For the attention of the patient

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during a medical examination.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

## 7 - FURTHER INFORMATION

The summary of safety and clinical performance of device shall be made available to the public via Eudamed (<https://ec.europa.eu/tools/eudamed>). Until the module is available on Eudamed, it can be accessed at [www.evolutisfrance.com](http://www.evolutisfrance.com).

For any further information please contact EVOLUTIS or its local representative.

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