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**GB - ORTHOPAEDIC IMPLANTS FOR SHOULDER ARTICULATION - UNIC STEMLESS**  
**INSTRUCTIONS FOR USE**  
**VERSION 2023-10**

For the attention of the surgeon and operation theatre staff

## **1. GENERAL**

The operating surgeon needs to be trained and familiar with the surgical technique, and must have read and understood the IFU (Instructions For Use) and the surgical technique provided by EVOLUTIS with the shipment of the material to the hospital, or consulted on-line or downloaded from the [www.evolutis-group.com](http://www.evolutis-group.com) site once the surgeon has requested and received a personal login code. Equally to the surgical procedure and to the positioning of the implants, the selection of the patient is important. An overweight or unsuitable functional requirement can create exceptional constraints and limit the lifespan of the arthroplasty. The warnings and the instructions must be strictly observed. Precautions and instructions should be rigorously followed. It is specified that all implants marketed by EVOLUTIS are supplied sterile (sterilization by gamma radiation).

## **2. INDICATIONS - CONTRA-INDICATIONS - COMPLICATIONS**

UNIC Stemless is indicated for the treatment of symptomatic pain and/or functional problems of the shoulder in patients whose skeleton is mature and only when pain killers and conservative treatment used correctly have failed. The anatomy and structure of the patient's articulation will need to be adapted to receive the selected implant(s).

### **2.1.Indications**

UNIC Stemless range is indicated for hemi or total anatomic shoulder arthroplasty.

The indications for UNIC Stemless are:

- Degenerative non-inflammatory arthropathy in shoulders presenting an intact or repairable rotator cuff: arthritis, avascular necrosis of the humeral head, post-traumatic arthritis.
- Inflammatory arthropathy such as rheumatoid arthritis.

Regardless of the pathology, the presence of a functional or repairable rotator cuff is absolutely necessary for the indication of a hemi or total anatomic shoulder arthroplasty.

### **2.2.Contra-indications**

Absolute contraindications to UNIC Stemless are: local or systemic infection, mental deficiency, neuromuscular disease, neurological or vascular disease, patients addicted to alcohol or psychotropic drugs and patients with allergies to prosthetic materials. The presence of a deficient or irreparable rotator cuff is an absolute contraindication to hemi or total anatomic shoulder arthroplasty.

Relative contra indications to UNIC Stemless are: excessive functional demands (sport with risk of falls, or with excessive functional demands beyond the limits of the prosthesis' resistance), overweight, insufficient bone stock or demineralization (osteoporosis) of the bone which would compromise the primary stability and the secondary fixation of the implant, drug abuse, pre-existing peri articular oncological pathology.

Implantation of prosthesis is not recommended in pregnant women.

### **2.3.Complications**

Implanting a prosthesis can lead to undesirable effects such as: infection, aseptic loosening, instability/dislocation, lysis/resorption, secondary deficiency of the rotator cuff, peri-prosthetic fracture, implant breakage/dissociation, glenoid erosion, pain, wear, progression of disease, arthrofibrosis, metal related pathology, stiffness, soft-tissue problems, osseous tissue problems, migration of the implant, malalignment/malposition/incorrect sizing of the implant, glenoid perforation, ossification, hematoma, wound problems and limb length discrepancy.

### **2.4.Benefits and performances**

Clinical benefits to the patients are pain relief, mobility restoration and function restoration.

Clinical performance and safety of the device are long term survival of the implant in situ with preservation of the function over time, and low frequency of complications.

### **3. PRECAUTIONS**

#### **3.1.Pre operative phase**

Prior to the operation, the surgeon must interview the patient about his physical and mental status and difficulties and inform him about the product and the different aspects of the operation. During the consultation the surgeon should explain the limits of the operation as well as the mechanical or the material limitations of the chosen implant(s). Factors which could limit the performance of the implant and alter its stability (eg. Overweight and functional demands) should be explained in order to reduce and delay the onset of complications linked to these factors. The necessity of following the surgeons' instructions and prescription must be explained to the patient. At the time of operation, the appropriate implants must be available.

#### **3.2.Manipulation**

For sterile implants the integrity of the packaging must be verified as it is a guarantee of sterility, as well as end of shelf-life date. If the packaging is damaged do not use, do not re-sterilize. Implants which have been unpacked and not implanted should not be re-sterilized. Implants must be manipulated with caution in order to avoid damage or scratching. They should be conserved in their original sealed packaging between 5 and 40°C. If protective components are delivered with the implant, they must not be removed until implantation of the implant. Implants are designed for single use and must never be re-implanted.

#### **3.3.Operative (Surgical) technique**

The surgeon must use known and recognized surgical techniques. Careful pre-operative planning on x rays should be done systematically. Templates are available for the majority of our implants.

#### **3.4.Fixation of implants**

Coated implants must not be implanted with cement. The appropriate implant of appropriate size should be used with the appropriate instrumentation and trials (if the latter are supplied).

#### **3.5.Post operative follow up**

Towards the end of the hospital stay patients should receive written advice and warnings about complimentary treatment, physio-therapy exercises and activity restrictions if necessary. In all cases validated post-operative follow up techniques must be used. The patient must be instructed to systematically inform that he or she has had an operation in case of subsequent MRI or scanner examination.

In the event of a serious incident, the user must notify the manufacturer and the competent authority.

#### **3.6.Removal of implant(s)**

In case of implant removal, the surgeon must use proven surgical techniques and appropriate instrumentation.

### **4. WARNINGS**

#### **4.1.Re-use of an implant**

Implants are single use and their re-use is forbidden. In case of re-use, the implant stability and its wear can be altered and the patient security compromised. If an implant is unwrapped by mistake or removed, they must be disposed of according to the legislation.

#### **4.2.Imaging exam of a patient with an Evolutis prosthesis**

Non-clinical testing has demonstrated that this prosthesis is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3T or less
- Maximum spatial gradient magnetic field of 19 T/m
- Maximum MR System reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the product is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

The appearance of an artifact caused by this device remains possible.

That's why the patient with a prosthesis must inform the healthcare staff before any MRI examination of the presence and site of implantation of his prosthesis.

### 4.3. Compatibility of implants and restrictions for use

Implants and trials provided by Evolutis should not be used with components supplied by another manufacturer since they might not be compatible.

### 4.4. Restrictions and specific recommendations for cobalt chrome (CoCr) components

The use of CoCr components should be avoided in:

- Patients with renal insufficiency
- Patients with cardiac insufficiency
- Patients of childbearing age.

For the other patients who will be implanted with a CoCr component, in case of suspected or proven failure of their prosthesis and whatever the component of the prosthesis involved in this failure, a medical follow-up including a cobalt blood test is recommended, knowing that cobalt is considered a CMR substance.

## 5. MATERIALS

The different materials used in UNIC Stemless implants are:

- Cobalt-chromium CoCr: ISO 5832-12 standard

ISO	C	Cr	Ni	Mo	Mn	N	Si	Fe	Balance
5832-12	0.15 to 0.35	26 to 30	≤ 1.0	5 to 7	≤ 1.0	≤ 0.25	≤ 1.0	≤ 0.75	Co

- Titanium alloy Ti6Al4V: ASTM F2924-14 standard

ASTM	C	N	Al	Y	O	V	Fe	H	Balance
F2924-14	≤ 0.08	≤ 0.05	5.50 to 6.75	≤ 0.005	≤ 0.20	3.50 to 4.50	≤ 0.30	≤ 0.015	Ti

- Calcium hydroxyapatite (coating): ISO 13779-2, 3 and 4 standards

- Ca10 (PO4)6 (OH)2 in compliance to ISO 13779-2, -3 and -4.

- Crystallinity: 45%

- Ca/P ratio: between 1.63 and 1.824

- CaO: less than 5%.

The compositions are given in mass fraction.

The material used for each individual device is specified on the packaging of each device.

The device:

- does not incorporate a substance which, if used separately, may be considered to be a medicinal product, including a medicinal product derived from human blood or human plasma.
- is not manufactured utilizing tissues or cells of human or animal origin, or their derivatives.
- is not composed of substances or combinations of substances that are intended to be introduced into the human body and that are absorbed by or locally dispersed in the human body.
- does not contain nanomaterials.

## 6. IMPLANT CARD

### 6.1. For the attention of the medical personnel

An implant card is provided in the prosthesis package. You should write the patient's name and the date of the operation and attach the patient label to the space provided. This card must be given to the patient along with the information specified below.

### 6.2. For the attention of the patient

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during an examination such as MRI, CT scan or ultrasound.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

## 7. FURTHER INFORMATION

The summary of safety and clinical performance of device shall be made available to the public via Eudamed (<https://ec.europa.eu/tools/eudamed>). Pending the availability of the module on Eudamed, it is available on the website [www.evolutis-group.com](http://www.evolutis-group.com).

For any further information please contact Evolutis or its local representative.

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